

Technical Data Sheet

Date: 11/07/2013

Product Code: LM-T1702

Description: Lymphosep, Lymphocyte Separation Media

FOR IN VITRO LABORATORY USE ONLY, NOT FOR DRUG, HUMAN VETERINARY USE

| | |
|----------------------------|---|
| Theoretical pH: | 7.0 ± 0.5 |
| Colour: | Colourless, clear solution |
| Storage conditions: | Room temperature |
| Shelf life: | 24 months |
| Sterility tests: | Bacteria in aerobic and anaerobic conditions Fungi and yeast |
| Density: | 1.077 ± 0.001 |
| Endotoxin: | < 10 EU/ml |

Recommended use:

Use aseptic technique when handling this product.

For in vitro laboratory use only, not for drug, human or veterinary use.

Application:

Lymphosep is designed for the simple, rapid isolation of lymphocytes from whole blood that has been diluted and treated with anti-coagulant or defibrinating agent.

For best results use blood drawn less than 2 hours before, do not use blood more than 24 hours from when it was drawn.

Uses:

- 1) Thoroughly mix the Lymphosep by inverting the bottle gently.
- 2) Aseptically transfer 3ml of Lymphosep to a 15ml centrifuge tube
- 3) Mix 2ml of defibrinated or heparined blood with 2ml of physiological saline (PBS w/o ca w/o Mg) or balanced salt solution (LM-S2041)
- 4) Carefully layer the diluted blood over 3ml of Lymphosep (room temperature) in a 15ml centrifuge, creating a sharp blood-Lymphosep interphase. DO NOT MIX! The quality of the separation is dependant upon a sharp interphase between the lymphocytes and the solution.
- 5) Centrifuge the tube at 400g at room temperature for 15 to 30 minutes. Centrifugation should sediment erythrocytes and polynuclear leukocytes and band mononuclear lymphocytes above the Lymphosep.
- 6) Aspirate the top layer of clear plasma to within the 2-3mm above the lymphocyte layer.
- 7) Aspirate the lymphocyte layer plus about half the Lymphosep layer below it and transfer it to a centrifuge tube. Add an equal volume of buffered balanced salt solution to the lymphocyte later in the centrifuge tube and centrifuge for 10 minutes at room temperature (18-25°C) at a speed sufficient to sediment to cells without damage i.e. 160-260g. Washing the cells removes Lymphosep and reduces the percentage of platelets.
- 8) Wash the cells again with a buffered balanced sale solution (LM-S2041) and re-suspend in the appropriate medium for your application.