

**Lymphosep, Lymphocyte Separation Media**

Code Produit : LM-T1702

<b><u>Theoretical pH</u></b> :	7.0 ± 0.5
<b><u>Osmolality</u></b> :	300 mOsm/kg ± 20
<b><u>Density</u></b> :	1.077 ± 0.001
<b><u>Colour</u></b> :	colourless, clear solution
<b><u>Storage conditions</u></b> :	Room temperature
<b><u>Shelf life</u></b> :	24 months
<b><u>Sterility tests</u></b> :	- bacteria in aerobic and anaerobic conditions - fungi and yeast
<b><u>Endotoxin</u></b> :	< 10 EU/ml
<b><u>Composition</u></b> :	Displayed on website and also available on request
<b><u>Recommended use</u></b> :	

- Respect storage conditions of the product
- Do not use the product after its expiry date
- Store product in an area protected from light (not necessary for saline solutions).
- Manipulate the product in aseptic conditions (e.g. : under laminar air flow)
- Wear clothes adapted to the manipulation of the product to avoid contamination (e.g. : gloves, mask, hygiene cap, overall...)

The product is intended to be used in vitro, in laboratory only. Do not use it in therapy, human or veterinary

applications.

## **Applications :**

Lymphosep is designed for the simple, rapid isolation of lymphocytes from whole blood that has been diluted and treated with anti-coagulant or defibrinating agent.

For best results, use blood drawn less than 2 hours before. Do not use blood more than 24 hours from when it was drawn.

## **Uses :**

- 1) Thoroughly mix the Lymphosep by inverting the bottle gently.
- 2) Aseptically transfer 3 ml of Lymphosep to a 15 ml centrifuge tube.
- 3) Mix 2 ml of defibrinated or heparinated blood with 2 ml of physiological saline (PBS w/o Ca w/o Mg) or balanced salt solution (L0615).
- 4) Carefully layer the diluted blood over 3 ml of Lymphosep (room temperature) in a 15 ml centrifuge, creating a sharp blood-Lymphosep interphase. DO NOT MIX! The quality of the separation is dependent upon a sharp interphase between the lymphocytes and the solution.
- 5) Centrifuge the tube at 400G at room temperature for 15 to 30 minutes. Centrifugation should sediment erythrocytes and polynuclear leukocytes and band mononuclear lymphocytes above the Lymphosep.
- 6) Aspirate the top layer of clear plasma to within 2-3 mm above the lymphocyte layer.
- 7) Aspirate the lymphocyte layer plus about half of the Lymphosep layer below it and transfer it to a centrifuge tube. Add an equal volume of buffered balanced salt solution to the lymphocyte layer in the centrifuge tube and centrifuge for 10 minutes at room temperature (18°C to 25°C) at a speed sufficient to sediment the cells without damage i.e., 160-260 g. Washing the cells removes Lymphosep and reduces the percentage of platelets.
- 8) Wash the cells again with buffered balanced salt solution (L0615) and resuspend in the appropriate medium for your applications.