

Versene

Product code:
LM-S2048

<u>Theoretical pH :</u>	7.2 ± 0.3
<u>Osmolality :</u>	280 mOsm/kg ± 10 %
<u>Colour :</u>	colourless, clear solution
<u>Storage conditions :</u>	+2°C to +8°C
<u>Shelf life :</u>	24 months
<u>Sterility tests :</u>	<ul style="list-style-type: none"> - Bacteria in aerobic and anaerobic conditions - Fungi and yeasts
<u>Endotoxin :</u>	< 1 EU/ml
<u>Activity test :</u>	Cells detachment test with the L929 cell line
<u>Composition :</u>	Available on request.
<u>Recommended use :</u>	<ul style="list-style-type: none"> - Respect storage conditions of the product - Do not use the product after its expiry date - Store product in an area protected from light (not necessary for saline solutions). - Manipulate the product in aseptic conditions (e.g. : under laminar air flow) - Wear clothes adapted to the manipulation of the product to avoid contamination (e.g. : gloves, mask, hygiene cap, overall...)

The product is intended to be used in vitro, in laboratory only. Do not use it in therapy, human or veterinary applications.

Description :

The Versene is used instead of Trypsin. It is a chelating agent that disperses the cells by cutting the cytoplasmic bridges between them.

The Versene has the advantage to bear the temperature of the autoclave, which is a guarantee of sterility. In addition, it is less aggressive cells than Trypsin, which is useful for studies on cell growth.

Uses :

To act, the Versene must be in a medium free of calcium and magnesium ions. It is generally used in PBS amended by deleting the calcium chloride and magnesium chloride.

- 1) Take 75 ml of Versene and complete to 100 ml with PBS. Autoclave at 120°C for 20 minutes.
- 2) Prepare a solution of CaCl₂ at 3.5 g/l and autoclave at 120°C for 20 minutes.
- 3) Discard the culture medium from the vial to transplant.
- 4) Rinse quickly but carefully the bottle with a little Versene-PBS, to remove calcium and magnesium ions remaining.
- 5) Add the Versene-PBS at 5 ml per 250 ml bottle.
- 6) Turn the bottle flat in the oven for 15 minutes.
- 7) After 15 minutes, take the bottle and add the CaCl₂ solution (0.2 ml per 250 ml) to neutralize the Versene.
- 8) Add the fresh medium directly under the same conditions than with Trypsin.
- 9) Shake, cells disperse.
- 10) Divide and put new bottles in the oven.

Signs of Deterioration :

Medium should be clear and free of particulate and flocculent material. Do not use if medium is cloudy or contains precipitate.