

---

**DMEM High Glucose w/ L-Glutamine w/o Sodium Pyruvate**

---

Product code :  
LM-D1109

<b><u>Theoretical pH:</u></b>	7.3 ± 0.3
<b><u>Osmolarity:</u></b>	334 mOsm/l ±10%
<b><u>Colour:</u></b>	Red solution
<b><u>Storage conditions:</u></b>	2 to 8° C
<b><u>Shelf life:</u></b>	12 months
<b><u>Endotoxin:</u></b>	<1 EU/ml
<b><u>Composition:</u></b>	Available on request
<b><u>Sterility Tests</u></b>	- Bacteria in aerobic and anaerobic conditions - Fungi and yeasts

**Cell Growth Test**

Medium tested for the ability to support cell growth with L929 cell line

**Recommended Use:**

- Respect storage conditions of the product
- Do not use the product after its expiry date
- Store product in an area protected from light (not necessary for saline solutions)
- Manipulate the product in aseptic conditions (e.g: under laminar air flow)
- Wear clothes adapted to the manipulation of the product to avoid contamination (e.g: gloves, mask, hygiene cap, overall...)

Product is provided for laboratory use only, and not for drug, human and veterinary use.

## **Description**

Lots of modifications of Eagle's medium have been developed since the creation of the first formulation. The most used Eagle's medium is the Dulbecco's Modified Eagle's Medium (DMEM).

It is a modification of Basal Medium Eagle (BME) that contains a concentration more important of amino acids and vitamins and also supplementary components. The original formulation contained 1000 mg/l of glucose and was used to culture embryonic mouse cells. The used of 4500 mg/l of glucose in the medium show an optimal cell growth for some cell lines.

## **Uses**

Supplements such as antibiotics should be added as sterile supplements to the medium. Storage conditions and shelf life of supplemented product will be affected by the nature of the supplements. Add 10ml/l of Sodium Pyruvate 100mM before using this medium.

## **Indication of Deterioration:**

Medium should be clear and free of particulate and flocculent material.

Do not use if the medium is cloudy or contains precipitate.

Other evidence of deterioration may include colour change or degradation of physical or performance characteristics.