

Technical Data Sheet

12/07/2017

RPMI 1640 Dutch Modification w/o L-Glutamine w/ 1g/l Sodium Bicarbonate w/20mM Hepes

Product Code :
LM – R1636

<u>Theoretical pH:</u>	7.3 ± 0.3
<u>Osmolality:</u>	300 mOsm/kg ± 10 %
<u>Colour:</u>	Clear orange solution
<u>Storage conditions:</u>	+2 to +8°C, protected from light
<u>Shelf life:</u>	24 months
<u>Sterility Tests:</u>	Bacteria aerobic-anaerobic Fungi/Yeast
<u>Endotoxin:</u>	<1 EU/ml
<u>Cell growth test:</u>	Medium tested for the ability to support SP2/0-Ag14 or Hela cell growth.
<u>Composition:</u>	Available on request

Recommended use:

- Respect storage conditions of the product
- Do not use the product after its expiry date
- Store product in an area protected from light (not necessary for saline solutions).
- Manipulate the product in aseptic conditions (e.g. : under laminar air flow)
- Wear clothes adapted to the manipulation of the product to avoid contamination (e.g. : gloves, mask, hygiene cap, overall...)

The product is intended to be used in vitro, in laboratory only. Do not use it in therapy, human or veterinary applications.

Applications:

RPMI 1640 medium has a broad spectrum of mammalian and hybridoma cell applications. It was developed by Moore and his co-workers at Roswell Park Memorial Institute in 1966 for the growth of human leukemia cells in monolayer or suspension cultures. It is typically supplemented with serum or serum substitutes.

Uses:

Supplements, such as antibiotics, should be added as sterile supplements to the medium. Storage conditions and shelf-life of supplemented product will be affected by the nature of the supplements.

Add 10.25ml/l of L-Glutamine 200mM before using this medium.

Signs of Deterioration:

Medium should be clear and free of particulate and flocculent material. Do not use if medium is cloudy or contains precipitate.

Other evidence of deterioration may include colour change or degradation of physical or performance characteristics.