
Phosphate Buffered Saline (PBS) w/o Calcium w/o Magnesium w/o Potassium Chloride

Product code:
LM-S2043

<u>Colour:</u>	colourless, clear solution
<u>Theoretical pH:</u>	7.4 ± 0.3
<u>Osmolality:</u>	290 mOsm/kg ± 10%
<u>Storage conditions:</u>	Room Temperature
<u>Shelf life:</u>	48 months
<u>Sterility tests:</u>	- Bacteria in aerobic and anaerobic conditions - Fungi and yeasts
<u>Endotoxin:</u>	< 1 EU/ml
<u>Composition:</u>	Displayed on website; also available on request

Recommended use:

- Respect storage conditions of the product
- Do not use the product after its expiry date
- Store product in an area protected from light (not necessary for saline solutions).
- Manipulate the product in aseptic conditions (e.g.: under laminar air flow)
- Wear clothes adapted to the manipulation of the product to avoid contamination (e.g.: gloves, mask, hygiene cap, overall...)

The product is intended to be used in vitro, in laboratory only. Do not use it in therapy, human or veterinary applications.

Application:

Phosphate-buffered saline (PBS) is a buffer solution used in biological research. It is a water-based salt solution containing sodium phosphate, sodium chloride and potassium phosphate. The osmolality and ion concentrations of the solutions are non-toxic to most cells.

Utilisation:

Supplements, such as antibiotics, should be added as sterile supplements to the buffer solution. Storage conditions and shelf-life of the supplemented product will be affected by the nature of the supplements.

Indications of deterioration:

Buffer solution should be clear and free of particulate and flocculent material.

Do not use if buffer solution is cloudy or contains precipitate.

Other evidence of deterioration may include degradation of physical or performance characteristics.